

INSTRUCTIONS QUICK START GUIDE

STEP 1

STEP 2

Fully Charge the Batteries before first use. (Approximately 6 hours, or when all 4 lights are on.)



Connect the battery pack to the Revitalizer cap. Be sure to locate and align the arrows.



STEP 3

Before use, disconnect the battery pack from the wall charger. If the battery pack is connected to the wall charger the device will not power on.



STEP 4

Place the Revitalizer cap on your head. THE REVITALIZER CAP WILL NOT FUNCTION UNLESS IT IS PLACED ON YOUR HEAD.

STEP 5

Locate the controller pod which is located on the cable between the Revitalizer cap and the battery pack. Press and hold the power button for 4 seconds. The cap is on when the white arch is illuminated. The laser treatment is ongoing when the yellow indicator light is illuminated.



Power is on when white arch is illuminated

STEP 6

To verify the Revitalizer Cap treatment has started, look into a mirror and you will notice that red laser lights pulsating on and off. The 30-minute treatment has started.

At the end of the 30-minute treatment cycle, the cap will automatically turn off.